

titramax VT

SULPHUR

Hydrogen sulphide and mercaptan sulphur in hydrocarbons and oils

Product description

The **Titramax VT SULPHUR** is designed for the determination of sulphur compounds as hydrogen sulphide and mercaptan sulphur in liquid hydrocarbons and oils.

The measurement uses a volumetric titration method with silver nitrate solution in an anhydrous medium. This measurement generates silver sulphide (Ag_2S) and silver mercaptides.

The device is conform to standards **ASTM D 3227, ISO 3012, UOP 163, UOP 212**.

The analysis must be done in an inert gas atmosphere to avoid an oxidation of sulphur compounds by air and thus false measurement results.

The titration speed is precisely adjusted to the reaction rate with control algorithms. Two end points get indicated – first for hydrogen sulphide and second for mercaptan sulphur.

At the end of the measurement, results are shown in mg/kg (ppm) hydrogen sulphide and/or mercaptan sulphur or several other units.



Titramax VT SULPHUR


Applications

The titrator is suitable for analysis of

- fuels (gasoline, kerosines, aviation turbine fuels)
- volatile distillates, middle distillates
- liquified petroleum gas (LPG)
- hydrocarbon gases
- liquid hydrocarbons
- naphtha (petroleum)
- low boiling oils
- light cycle oils
- heavy cycle oils

Advantages

- Complete measuring system for the determination of hydrogen sulphide and mercaptan sulphur
- Fully-automatic volumetric titration
- Precise adjustment of the titration parameters by control algorithms
- Preset measurement method allows an immediate start
- The result output can be adjusted to your needs by using a formula generator



Titration graph of a sample

Features

The **Titramax VT SULPHUR** consists of


- an automatic volumetric titrator with potentiometric indication
- a titration vessel with stirrer unit

The determination of hydrogen sulphide/mercaptan sulphur is based on

- a potentiometric titration with silver nitrate solution in an anhydrous medium
- a precise indication by a selective silver electrode, which is stable over long periods

Steps of the analysis are

1. Standardization of the titration solution
2. Titration of oil sample or absorption liquid



Technical specifications

Measurement method:

Types of result:

Measuring range:

Display resolution:

Measurement range pH / mV:

Display resolution pH / mV:

Accuracy pH / mV (without sensor):

Measurement range μ A:

Display resolution μ A:

Accuracy μ A (without sensor):

Measurement range temperature $^{\circ}$ C:

Amplifier input impedance:

Burette resolution:

Dosing accuracy according DIN EN ISO 8655, part 3: Accuracy 0.15 % / Precision 0.05 - 0.07 % (depending on the used exchange unit)

Filling time:

Power supply:

Power input:

Stirrer connection:

Dimensions:

Weight:

Volumetric titration

mg/kg (ppm) H_2S and/or mercaptan sulphur, formula generator available

Hydrogen sulphide/sulphide: less than 1 to several thousand mass-ppm

Mercaptan sulphur: 3 - 500 mass-ppm typically

0.01 ppm

- 3.0 ... 18.00 / - 2000 ... 2000

0.001 / 0.1

0.002 / 0.1 mV \pm 1 digit

0 ... 100

0.1

0.2 \pm 1 digit

- 75 ... 175

> 1 $\cdot 10^{13}$ ohms

10,000 steps for 10 mL / 20 mL \pm 0.15 %

20 sec

External plug-in power supply 100 - 240 V, 50/60 Hz

30 VA

12 V DC out, 500 mA

30 x 45 x 30 cm (W x H x D), height with exchange unit

Approx. 3.5 kg (with exchange unit and empty reagent bottle)

ECH Elektrochemie Halle GmbH

Otto-Eißfeldt-Str. 8
D-06120 Halle (Saale)
Germany

Tel.: +49 (0) 345 279570-0

Fax: +49 (0) 345 279570-99

Email: info@echscientific.com • www.ech.de • www.aquamaxkf.com

ECH Scientific Limited

Building 69, Wrest Park, Silsoe
Bedfordshire, MK45 4HS
United Kingdom

Tel.: +44 (0) 1525 404747

Fax: +44 (0) 1525 404848



SCIENTIFIC

part of ECH Elektrochemie Halle

Global Sales Division

the ECH advantage

in-lab | mobile | on-line | process