

# Manual Oil Sampling Kit

## Product description

During the sampling of oils you must avoid any contact with ambient air.

The Manual Oil Sampling Kit is designed for contamination-free taking of oil samples.

A constant temperature of oil sample is important for correct analysis. The temperature is automatically measured during flow of oil through the Manual Oil Sampling Kit. When the temperature of oil sample is constant you can fill the sample into the gas tight glass syringe.

The Manual Oil Sampling Kit is easy to use. Versions for various syringe sizes are available (10 mL, 50 mL, 100mL).

## Applications

- Mineral oils
- Brake fluids
- Biological oils
- Food- and pharmaceutical products
- Vegetable oils
- Oxidation sensitive oils

## Advantages

- Gas tight sampling kit for oils
- No contact of oil sample with ambient air
- Versions for various syringe sizes (10 mL, 50 mL, 100 mL)



Manual Oil Sampling Kit with glass syringe (100 mL)

## We are here for you



ECH Elektrochemie Halle GmbH  
Otto-Eissfeldt-Str. 8  
D-06120 Halle (Saale)  
Germany  
Tel.: +49 345 279570-0  
Fax: +49 345 279570-99  
E-mail: [info@ech.de](mailto:info@ech.de)  
Website: [www.ech.de](http://www.ech.de)